

# Tiras con esporas Spore Strips Tiras com esporos

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# Spore Strips

BT40

For Ethylene Oxide and Dry Heat sterilization

## Composition

Each Bionova® BT40 packet consists of a paper strip soaked with *Bacillus atrophaeus* ATCC 9372 spore population.

## Product Description

The packet containing Bionova® BT40 spore strip is specifically designed to control either ethylene oxide (EO) or dry heat (Dry) sterilization processes. Spore strips should be incubated in Bionova® MC1030 growth medium tube or in growth medium appropriate for *Bacillus atrophaeus* growth.

If sterilization process was not successful, MC1030 growth medium will turn to yellow after incubation at 37±2°C, thus indicating the existence of living *Bacillus atrophaeus* spores in spore strips.

If sterilization process was correct, MC1030 indicating medium will keep its original color. A last reading should therefore be done after 48 hours of growth medium incubation at 37±2°C.

## Warning!

Do not use BT40 spore strip packets to control steam sterilization cycles, chemical vapors, radiation or other sterilization processes. Do not re-use spore strip packets.

## Important Considerations

It is advisable to use at least two BT40 spore strip packets to control each sterilization cycle.

## Storage

Store in a dark place and at temperatures between 15-30°C, 35-60 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

## Directions for use

1. Identify Bionova® BT40 packet by writing on its back the sterilizer number (in case it has more than one), load number and processing date.

2. Pack it in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the packet in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.

3. Sterilize as usual.

4. After sterilization process has finished you can:

A. Remove Bionova® BT40 packet from the sterilization package used for processing and incubation.

B. For EO sterilization processes you can ventilate package containing Bionova® BT40 packet along with the remaining load or you can remove it before ventilation process for spore strips incubation in MC1030 growth medium or any other appropriate growth medium.

5. Tear up the packet at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer spore strip with a sterile clamp. **Important:** Use latex gloves and chin strap when transferring spore strips from the packet to MC1030 growth medium tube. **Warning!** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube and of Bionova® BT40 packet.

6. Incubate spore strips at 37±2°C.

**Important:** Use a spore strip which has not been under sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

7. Incubate processed spore strips along with the strip used as positive control during a maximum of 48 hours at 37±2°C. It is advisable to make observations every 10 hours. Growth medium indicator color turning to yellow shows a failure on sterilization process. If after 48 hours no color change is visible on MC1030 growth medium, the result is negative (i.e. the sterilization process was successful). MC1030 indicating medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discharge them as follows: **Warning!** Do not re-use sterilizer until spore strip growth result turns negative (MC1030 growth medium containing processed spore strips remains in its original color).

## Disposal

Dispose of the growth medium tubes already used to analyse growth of the spores inside Bionova® BT40 strips according to health regulations in your country. Positive MC1030 growth medium tubes can be sterilized in steam sterilizer at 121°C at least during 15 minutes, or at 132°C for 10 minutes in gravity displacement steam sterilizer, or at 132°C for 4 minutes in vacuum steam sterilizer.

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