

Tiras con esporas Spore Strips Tiras com esporos

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
Parameters determined at time of manufacture according to ISO 11138 (Parts 1 and 3) and IRAM 37102 (Parts 1 and 3) standards. The showed values are reproducible only under the same conditions under which they were determined.

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

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Industria Argentina
Made in Argentina

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For Steam sterilization

Composition

Each Bionova® BT50 packet consists of a paper strip soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product Description

The packet containing Bionova® BT50 spore strip is specifically designed to control steam (Steam) sterilization processes. Spore strips should be incubated in Bionova® MC20 growth medium tube or in growth medium appropriate for *Geobacillus stearothermophilus* growth.

If sterilization process was not successful, MC20 growth medium will turn to yellow after incubation at 56±2°C, thus indicating the existence of living *Geobacillus stearothermophilus* spores in spore strips.

If sterilization process was correct, MC20 indicating medium will keep its original color. A last reading should therefore be done after 24 hours of growth medium incubation at 56±2°C.

Warning!

Do not use BT50 spore strip packets to control EO or dry heat sterilization cycles, radiation or other sterilization processes. Do not re-use spore strip packets.

Important Considerations

It is advisable to use at least two BT50 spore strip packets to control each sterilization cycle.

Storage

Store in a dark place and at temperatures between 15-30°C, 35-60 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Identify Bionova® BT50 packet by writing on its back the sterilizer number (in case it has more than one), load number and processing date.

2. Pack it in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the packet in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.

3. Sterilize as usual.

4. After sterilization process has finished, remove Bionova® BT50 packet from the package used for processing and incubation.

5. Tear up the packet at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer spore strip with a sterile clamp to MC20 growth medium tube or to any other appropriate growth medium. **Important!** Use latex gloves and chin strap when transferring spore strips from the packet to growth medium tube. **Warning!** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube and of Bionova® BT50 packet.

6. Incubate spore strip at 56±2°C.

Important! Use a spore strip which has not been under sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

7. Incubate processed spore strips along with the strip used as positive control during a maximum of 24 hours at 56±2°C. It is advisable to make observations every 10 hours. Growth medium indicator color turning to yellow shows a failure on sterilization process. If after 24 hours no color change is visible on MC20 growth medium, the result is negative (i.e. the sterilization process was successful). MC20 indicating medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discharge them immediately as indicated below. **Warning!** Do not re-use sterilizer until spore strip growth result turns negative (MC20 growth medium containing processed spore strip remains in its original color).

Disposal

Dispose of the growth medium tubes already used to analyse Bionova® BT50 spore strips growth according to health regulations in your country. Positive MC20 growth medium tubes can be sterilized in steam sterilizer at 121°C at least during 15 minutes, or at 132°C for 10 minutes in gravity displacement steam sterilizer, or at 132°C for 4 minutes in vacuum steam sterilizer.

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