

# Sterilizing gas ampoules

## BIOLENE® Sterilization Kit

This method of sterilization is so simple that can be carried out anywhere. No gas or air installation, electrical connections or specially designed rooms are required, making it the best option for areas under emergency conditions.

The Sterilization Kit consists of one sterilizing gas ampoule (ethylene oxide + inert gases) wrapped in a diffusing membrane, one liner bag and one closure to tighten the bag.

Any kind of material, previously packed, can be effectively sterilized with BIOLENE® Sterilization Kit. Materials are not damaged, since this sterilization system does not use heat and/or pressure.



CODE	AMPOULE		BAG		PACKAGING	
	Volumen (cm <sup>3</sup> )	ETO (%)	Measures (cm)	Capacity (litres)	Standard (Kits)	Export (Kits)
BL-05	5	90	31 x 60	11	box x 60	box x 20
BL-20	20	97	56 x 92	34	box x 20	box x 5

Products manufactured in Argentina



Address:  
Constitución 2610  
(C1254AAR)  
Buenos Aires  
Argentina

Telephone:  
(+54 11) 4308.4963

Fax:  
(+54 11) 4941.0346

e-mail:  
info@biolene.com

www.biolene.com

Quality Management System



CERTIFIED according to ISO 9001:2000

**BIOLENE®**  
ETHYLENE OXIDE STERILIZATION