

Ethylene Oxide cartridges

Ethylene Oxide cartridges

- 100% ethylene oxide.
- Single-use cartridge: after the sterilization (and aeration) cycle is complete, it is discarded with non-incinerated waste.
- Highly effective against bacteria, fungus and virus.
- Compatible with a wide range of sterilizers available in the market.

CODE	MEASURES (mm)		CAPACITY	
	Diameter	Height	Grams	Ounces
BL-07	25	85	10	0,36
BL-15	25	85	15	0,53
BL-25	25	105	30	1,06
BL-50	35	105	50	1,77
BL-67	35	105	67	2,37
BX-100	38	165	100	3,52
BX-134	38	165	127	4,48
BX-170	45	180	170	5,99
BQ-70	35	140	70	2,47
BQ-120	50	133	120	4,24
BQ-130	50	133	130	4,59
BQ-200	50	174	200	7,06



Biolene has registered at EPA (U. S. Environmental Protection Agency) its metallic cartridges and establishment.



Products manufactured in Argentina

Address:
Constitución 2610
(C1254AAR)
Buenos Aires
Argentina

Telephone:
(+54 11) 4308.4963

Fax:
(+54 11) 4941.0346

e-mail:
info@biolene.com

www.biolene.com

Quality Management System



CERTIFIED according to ISO 9001:2000



BIOLENE®
ETHYLENE OXIDE STERILIZATION